

# Hardware Code

```
#include <Arduino.h>
#include <bluefruit.h>
#include <movingAvg.h>
#include <arduinoFFT.h>
#include <Adafruit_LSM6DS3.h> // For the onboard IMU

// BLE Configuration
BLEService healthService("1810");
BLECharacteristic hrChar("2A37", BLERead | BLENotify, 2);
BLECharacteristic spO2Char("2A5F", BLERead | BLENotify, 2);
BLECharacteristic eegChar("2A38", BLERead | BLENotify, 20); // EEG data
packet
BLECharacteristic battChar("2A19", BLERead | BLENotify, 1);

// Pins
#define VIBRATION_PIN D2
#define BATTERY_PIN A0
#define EEG_INPUT_PIN A1 // BioAmp EXG Pill output

// EEG Configuration
const uint16_t EEG_SAMPLES = 128; // Must be power of 2
const float SAMPLING_RATE = 250.0; // Hz
ArduinoFFT<double> FFT;
double vReal[EEG_SAMPLES];
double vImag[EEG_SAMPLES];
unsigned long samplingPeriod;

// Moving averages
movingAvg hrAvg(5);
movingAvg spO2Avg(5);
movingAvg eegAvg(5);

// Battery params
```

```

const float BAT_MAX = 3.7f;
const float BAT_MIN = 3.0f;

// EEG Band Definitions (Hz)
enum EEGBand {
    DELTA = 0, // 0.5-4Hz
    THETA,    // 4-8Hz
    ALPHA,    // 8-12Hz
    BETA,     // 12-30Hz
    GAMMA,    // 30-100Hz
    BAND_COUNT
};

enum EmergencyType {
    HR_EMERGENCY,
    SPO2_EMERGENCY,
    SPO2_CRITICAL,
    EEG_ARTIFACT_DETECTED,
    SEIZURE_PATTERN,
    FALL_DETECTED
};

// Struct for EEG analysis
struct EEGAnalysis {
    uint8_t bandPercents[BAND_COUNT];
    uint8_t dominantBand;
};

Adafruit_LSM6DS3 lsm6ds3; // IMU object
// Fall detection parameters
const float FALL_ACCEL_THRESHOLD = 2.5f; // g-force threshold (2.5g)
const float FREE_FALL_THRESHOLD = 0.5f; // g-force threshold for free-fall
//
const unsigned long IMPACT_DURATION = 1000; // ms to detect post-impact stillness
const float POST_FALL_ANGLE = 45.0f; // degrees from vertical

```

```

// Variables for fall detection
unsigned long fallStartTime = 0;
bool potentialFall = false;
float preFallOrientation[3] = {0};

void setup() {
  Serial.begin(115200);
  while(!Serial) delay(10);

  // Initialize hardware
  pinMode(VIBRATION_PIN, OUTPUT);
  pinMode(EEG_INPUT_PIN, INPUT);

  hrAvg.begin();
  spO2Avg.begin();
  eegAvg.begin();

  // Calculate sampling period
  samplingPeriod = 1000000 / SAMPLING_RATE;

  // Initialize BLE
  Bluefruit.begin();
  Bluefruit.setName("MindHealth_Tracker");
  Bluefruit.setTxPower(4);

  setupBLE();
  Serial.println("BLE Ready");
}

void setupBLE() {
  healthService.begin();

  hrChar.setProperties(BLERead | BLENotify);
  hrChar.setPermission(SECMODE_OPEN, SECMODE_NO_ACCESS);
  hrChar.begin();

  spO2Char.setProperties(BLERead | BLENotify);
  spO2Char.setPermission(SECMODE_OPEN, SECMODE_NO_ACCESS);

```

```

spO2Char.begin();

eegChar.setProperties(BLERead | BLENotify);
eegChar.setPermission(SECMODE_OPEN, SECMODE_NO_ACCESS);
eegChar.begin();

battChar.setProperties(BLERead | BLENotify);
battChar.setPermission(SECMODE_OPEN, SECMODE_NO_ACCESS);
battChar.begin();

Bluefruit.Advertising.addService(healthService);
Bluefruit.Advertising.start();
}

void loop() {
    static uint32_t lastSampleTime = 0;
    static uint16_t sampleCount = 0;

    // // BLE maintenance
    // Bluefruit.process();

    // EEG Sampling
    if (micros() - lastSampleTime >= samplingPeriod && sampleCount < EEG_
SAMPLES) {
        vReal[sampleCount] = analogRead(EEG_INPUT_PIN);
        vImag[sampleCount] = 0;
        sampleCount++;
        lastSampleTime = micros();
    }

    // Process EEG when buffer full
    if (sampleCount >= EEG_SAMPLES) {
        EEGAnalysis analysis = processEEG();
        sendEEGData(analysis);
        sampleCount = 0;
    }

    // Simulate HR/SpO2 (replace with actual MAX30102 code)

```

```

if (millis() % 1000 == 0) {
    updateHRSpO2();
    checkBattery();
}

checkEmergencies();
}

EEGAnalysis processEEG() {
    EEGAnalysis result = {{0}, 0};

    // Apply FFT
    FFT.windowing(vReal, EEG_SAMPLES, FFT_WIN_TYP_HAMMING, FFT_FORWARD);
    FFT.compute(vReal, vImag, EEG_SAMPLES, FFT_FORWARD);
    FFT.complexToMagnitude(vReal, vImag, EEG_SAMPLES);

    // Classify frequency bands
    float bandTotals[BAND_COUNT] = {0};
    float totalPower = 0;

    for (int i = 1; i < EEG_SAMPLES/2; i++) {
        float freq = i * (SAMPLING_RATE / EEG_SAMPLES);
        float magnitude = vReal[i];

        if (freq < 4) bandTotals[DELTA] += magnitude;
        else if (freq < 8) bandTotals[THETA] += magnitude;
        else if (freq < 12) bandTotals[ALPHA] += magnitude;
        else if (freq < 30) bandTotals[BETA] += magnitude;
        else if (freq < 100) bandTotals[GAMMA] += magnitude;

        totalPower += magnitude;
    }

    // Calculate percentages
    if (totalPower > 0) {
        for (int i = 0; i < BAND_COUNT; i++) {
            result.bandPercents[i] = (bandTotals[i] / totalPower) * 100;
        }
    }
}

```

```

    }

    // Find dominant band
    float maxVal = 0;
    for (int i = 0; i < BAND_COUNT; i++) {
        if (bandTotals[i] > maxVal) {
            maxVal = bandTotals[i];
            result.dominantBand = i;
        }
    }
}

return result;
}

void sendEEGData(EEGAnalysis analysis) {
    uint8_t eegData[20] = {0};

    // Format: [dominantBand, delta%, theta%, alpha%, beta%, gamma%]
    eegData[0] = analysis.dominantBand;
    for (int i = 0; i < BAND_COUNT; i++) {
        eegData[i+1] = analysis.bandPercents[i];
    }

    eegChar.write(eegData, sizeof(eegData));
}

void updateHRSpO2() {
    // Simulate data (replace with MAX30102 readings)
    uint16_t hr = 70 + random(-5, 5);
    uint16_t spo2 = 97 + random(-1, 1);

    hrAvg.reading(hr);
    spO2Avg.reading(spo2);

    uint16_t hrValue = hrAvg.getAvg();
    uint16_t spO2Value = spO2Avg.getAvg();
    hrChar.write(&hrValue, sizeof(uint16_t));
}

```

```

    spO2Char.write(&spO2Value, sizeof(uint16_t));

}

void checkBattery() {
    int raw = analogRead(BATTERY_PIN);
    float voltage = (raw * 3.3f / 1024.0f) * 2; // Voltage divider
    uint8_t level = map(voltage * 100, BAT_MIN * 100, BAT_MAX * 100, 0, 100);

    battChar.write(&level, 1);

    if (level < 20) {
        triggerVibration(2, 500); // Double buzz for low battery
    }
}

void checkEmergencies() {
    static uint32_t lastCheckTime = 0;
    if (millis() - lastCheckTime < 1000) return; // Run once per second
    lastCheckTime = millis();

    // Current vital readings
    float currentHR = hrAvg.getAvg();
    float currentSpO2 = spO2Avg.getAvg();

    // Threshold values (customize as needed)
    const float HR_CRITICAL_LOW = 40.0f;
    const float HR_CRITICAL_HIGH = 120.0f;
    const float SpO2_CRITICAL = 90.0f;
    const float EEG_ARTIFACT_THRESHOLD = 500.0f; // Raw EEG amplitude threshold

    // Check heart rate emergencies
    if (currentHR < HR_CRITICAL_LOW || currentHR > HR_CRITICAL_HIGH) {
        triggerEmergencyAlert(HR_EMERGENCY);
        Serial.print("HR Emergency: ");
        Serial.println(currentHR);
    }
}

```

```

// Check SpO2 emergencies
if (currentSpO2 < SpO2_CRITICAL) {
    triggerEmergencyAlert(SPO2_EMERGENCY);
    Serial.print("SpO2 Emergency: ");
    Serial.println(currentSpO2);
}

// EEG Spike Detection (simple amplitude check)
for (int i = 0; i < EEG_SAMPLES; i++) {
    if (abs(vReal[i]) > EEG_ARTIFACT_THRESHOLD) {
        triggerEmergencyAlert(EEG_ARTIFACT_DETECTED);
        Serial.println("EEG Artifact Detected");
        break;
    }
}

// Check for seizure-like patterns (high gamma dominance)
EEGAnalysis latestAnalysis = processEEG();
if (latestAnalysis.bandPercents[GAMMA] > 40 &&
    latestAnalysis.bandPercents[BETA] > 30) {
    triggerEmergencyAlert(SEIZURE_PATTERN);
    Serial.println("Potential seizure pattern detected");
}

// Fall detection using IMU (if available)
if (isFallingDetected()) {
    triggerEmergencyAlert(FALL_DETECTED);
    Serial.println("Fall detected");
}

void triggerEmergencyAlert(EmergencyType type) {
    switch(type) {
        case HR_EMERGENCY:
            // Triple vibration for heart emergency
            triggerVibration(3, 300);
            break;
    }
}

```



```

case SPO2_CRITICAL:
    // Double long vibration for SpO2
    triggerVibration(2, 500);
    break;

case EEG_ARTIFACT_DETECTED:
    // Quick bursts for EEG issues
    triggerVibration(4, 100);
    break;

case SEIZURE_PATTERN:
    // Continuous vibration for seizure
    digitalWrite(VIBRATION_PIN, HIGH);
    delay(3000); // 3 second continuous
    digitalWrite(VIBRATION_PIN, LOW);
    break;

case FALL_DETECTED:
    // Alternating pattern for fall
    for (int i = 0; i < 5; i++) {
        digitalWrite(VIBRATION_PIN, HIGH);
        delay(100);
        digitalWrite(VIBRATION_PIN, LOW);
        delay(100);
    }
    break;
}
}

void setupIMU() {
    if (!lsm6ds3.begin_I2C()) {
        Serial.println("Failed to initialize LSM6DS3 IMU");
        while (1);
    }

    // Configure IMU settings
    lsm6ds3.setAccelRange(LSM6DS3_ACCEL_RANGE_16_G);

```

```

lsm6ds3.setAccelDataRate(LSM6DS_RATE_104_HZ);
lsm6ds3.setGyroRange(LSM6DS_GYRO_RANGE_2000_DPS);
lsm6ds3.setGyroDataRate(LSM6DS_RATE_104_HZ);
}

bool isFallingDetected() {
    sensors_event_t accel;
    sensors_event_t gyro;
    sensors_event_t temp;
    lsm6ds3.getEvent(&accel, &gyro, &temp);
    // Calculate total acceleration magnitude
    float accelMag = sqrt(sq(accel.acceleration.x) +
                          sq(accel.acceleration.y) +
                          sq(accel.acceleration.z));
    // Detect free-fall (low g-force)
    if (accelMag < FREE_FALL_THRESHOLD && !potentialFall) {
        potentialFall = true;
        fallStartTime = millis();

        // Store pre-fall orientation
        preFallOrientation[0] = accel.acceleration.x;
        preFallOrientation[1] = accel.acceleration.y;

        // Check post-impact stillness
        if (checkPostFallStillness()) {
            // Verify orientation change
            float currentOrientation[3] = {
                accel.acceleration.x,
                accel.acceleration.y,
                accel.acceleration.z
            };

            float angleChange = calculateOrientationChange(preFallOrientation, currentOrientation);

            if (angleChange > POST_FALL_ANGLE) {
                potentialFall = false;
            }
        }
    }
}

```

```

        return true;
    }
}
}
// Reset if too much time has passed
if (potentialFall && (millis() - fallStartTime) > 2000) {
    potentialFall = false;
}
return false;
}

```

```

bool checkPostFallStillness() {
    unsigned long startTime = millis();
    float initialAccel[3];

    sensors_event_t accel;
    lsm6ds3.getEvent(&accel, NULL, NULL);

    initialAccel[1] = accel.acceleration.y;
    initialAccel[2] = accel.acceleration.z;

    while (millis() - startTime < IMPACT_DURATION) {
        lsm6ds3.getEvent(&accel, NULL, NULL);

        // Check for significant movement
        float diff = sqrt(sq(accel.acceleration.x - initialAccel[0]) +
                        sq(accel.acceleration.y - initialAccel[1]) +
                        sq(accel.acceleration.z - initialAccel[2]));

        if (diff > 0.3f) { // Movement threshold
            return false;
        }

        delay(50);
    }

    return true;
}

```

```

}

float calculateOrientationChange(float before[3], float after[3]) {
    // Calculate dot product
    float dot = before[0]*after[0] + before[1]*after[1] + before[2]*after[2];

    // Calculate magnitudes
    float magBefore = sqrt(sq(before[0]) + sq(before[1]) + sq(before[2]));
    float magAfter = sqrt(sq(after[0]) + sq(after[1]) + sq(after[2]));

    // Calculate angle in radians then convert to degrees
    float angle = acos(dot / (magBefore * magAfter)) * 57.2958;
    return angle;
}

void triggerVibration(uint8_t count, uint16_t duration) {
    for (uint8_t i = 0; i < count; i++) {
        digitalWrite(VIBRATION_PIN, HIGH);
        delay(duration);
        digitalWrite(VIBRATION_PIN, LOW);
        if (i < count-1) delay(200);
    }
}

```